# Flexible Endoscope Cleaning, Disinfecting, Storage and Tips


For detailed steps on cleaning and disinfecting your endoscope, refer to the IFU that was provided with your device.

**NOTE:** Use of proper personal protective equipment during endoscope handling/reprocessing is recommended.

**\*Solutions NOTE:** Follow manufacturer's instructions when preparing cleaning and disinfection solutions for soak times, temperature and rinsing. The supplied IFU lists approved method for high level disinfection.


## 1. Pre-cleaning

Pre-clean immediately after every patient procedure, this will prevent patient material from adhering to the endoscope. Dried material is difficult to remove and can render the disinfection or sterilization process ineffective. Gently wipe all debris from the insertion tube and distal tip using a soft, lint-free cloth and water with a low sudsing detergent\* making sure to also clean the lens at the distal tip face.

**CAUTION:** Use care when wiping down the insertion tube and cleaning the distal tip.


**CAUTION:** Do not use excessive force or a twisting motion when wiping down the insertion tube. Use firm but gentle pressure paying specific attention to the distal lens.

## 2. Leak Testing

After every procedure, in preparation for high-level disinfection or sterilization, the endoscope should first be tested to ensure that it is watertight by leak testing.

**DISCONTINUE USE OF THE ENDOSCOPE IF A LEAK IS DETECTED**

Connect the leak tester by aligning the pin in the Leak Tester connector with the groove on the Vent, pushing down, and rotating clockwise.

Pump the hand bulb of the leak tester until the needle is within the 140-180mmHg. Observe the needle for 10 seconds and for any drop in pressure. **Do not exceed 180mmHg**

Immerse the entire endoscope in water, observe for 30 seconds. A steady stream of bubbles indicates a leak.

Absence of air bubbles confirms that the endoscope is watertight, remove endoscope from water and disconnect the leak tester - Proceed to Step 3.

## 3. Cleaning

1. Fill a basin with a freshly made solution of a low sudsing detergent\*.
2. Immerse the endoscope. Wash all debris from the exterior by gently brushing and wiping.
3. Rinse the endoscope in clean water.
4. Wipe dry the exterior of the endoscope with a soft, clean cloth.

**NOTE:** Before immersing scope, ensure Vent Cap is OFF

## 4. Disinfecting

1. Prepare the high level disinfectant solution\*.
2. Completely immerse the endoscope in the high level disinfectant solution\*.
3. Thoroughly rinse the exterior of the endoscope with clean water\*.
4. Dry the exterior of the endoscope with a soft, clean towel.
5. Ensure the Eyepiece and Distal Tip lenses are free from solution residue.

**NOTE:** Before immersing scope, ensure Vent Cap is OFF

## 5. Storage

- The endoscope should be stored in a dry, well-ventilated environment.
- Do not store the endoscope in the carrying case. The carrying case is not sufficiently ventilated for storage of the endoscope. The carrying case is intended only for transport of the endoscope.
- The endoscope should be completely clean and dry before storing.
- When storing the endoscope, be sure to keep the insertion shaft straight

**NOTE:** Ensure Vent Cap is OFF  
**NOTE:** Recommended storage is the EndoCaddy

## Notes:

Routine Maintenance: Optim recommends replacing the O-ring with the provided lubricated o-ring every 6 months or 300 disinfectant cycles.

Optional Lubricant: Apply approved lubricant to the o-ring using your fingers.

Refer to IFU for detailed instructions.

Please contact Optim's Customer Service Department with any questions at (800) 225-7486.

64 Technology Park Road, Sturbridge, MA  
sales@optim-llc.com / www.optim-llc.com  
800.225.7486 / 1.508.347.5100

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