

Rigid Sinuscope & PLS Bundle



Our Rigid Sinuscope combined with the Portable LED Light Source (PLS) are designed and manufactured to the highest quality standards to ensure reliable performance to meet the rigors of real-world office and surgical use. The portable solution eliminates the need for external light sources and lamp replacement.

The rigid sinuscope combines precision-machined components with superior optics to provide exceptional clarity and brightness.

Features & Benefits

- “Universal” eyepiece compatible with virtually all add-on camera systems
- All Optim ENT Rigid Sinuscopes are soakable and autoclavable
- Utilization with the Optim PLS eliminates cables, providing freedom of movement during the exam
- The patented optical coupling technology of the Optim PLS focuses and shapes the light, delivering superior image brightness and optical clarity
- Rechargeable batteries provide power for an entire day of exams.



Bundle includes (1) Rigid Sinuscope, Portable LED Light Source with (2) batteries and charger all in a compact carrying case.

Specifications

Model	Diameter	Working Length	Viewing Direction
RPK 4000	4.0 mm	175 mm	0°
RPK 4030	4.0 mm	175 mm	30°
RPK 4070	4.0 mm	175 mm	70°
RPK 2700	2.7 mm	187 mm	0°
RPK 2730	2.7 mm	187 mm	30°
RPK 2770	2.7 mm	187 mm	70°



64 Technology Park Road, Sturbridge, MA 01566
Toll-free: 800.225.7486 Worldwide: 1.508.347.5100
Fax: 508.347.2380
Email: sales@optim-llc.com
Website: www.optim-llc.com



Optim is ISO 9001:2008 / ISO13485:2003 Certified