



High Performance Diagnostics in a Simple and Mobile System



Designed to optimize the freedom to perform the Flexible Endoscopic Evaluation of Swallowing (FEES) procedure in a variety of settings, the ENTity™ MobileFEES system enables procedures from patient bedside, to office settings, or throughout the hospital. The ENTity MobileFEES system is the optimal choice for nursing homes, rehabilitation, hospitals, clinics and independent Speech Language Pathologists.

Training & Service

- On-site equipment training
- Superior customer service
- After-install support

Portability

- Portable and lightweight case
- Bring the scope to the patient
- Use in the hospital / bedside / office

User Friendly Software

- Capture still images & streaming videos
- EMR compatible & HIPAA compliant
- Windows based

MobileFEES System Includes:

- ENTity XL or SDXL Laryngoscope with Internal LED Light Source
- 15.6 Dell Touchscreen Laptop - Windows 8 Professional
- VidCap USB Video Camera Kit with Sopro Software
- Pediatric and Adult FEES Report Templates for Documentation
- Foot Pedal Camera Activator
- 2TB (Terabyte) Portable External Hard Drive Storage
- HP Officejet 100 Color Mobile Printer
- Two-Year Manufacturers Warranty on ENTity XL / SDXL *with Free Loaner Equipment*
- Toll-Free Clinical Support Hot-line
- Software Support with remote access - Available: M-F / 9-7 EST
- On-Site Equipment Training by our Clinical Specialist
- Accessories: Battery Charger, (4) Rechargeable Lithium Ion Batteries, Leak Tester & Carrying Case



The MobileFEES System can utilize both the ENTity XL flexible laryngoscope, a standard 3.6mm for traditional anatomy, or the ENTity SDXL with a slim, small 2.7mm diameter for patients with smaller or more difficult anatomy

