

# Flexible Nasopharyngoscope Model ENT 3.6F

*Precision endoscopy with proven performance in two flexible configurations*



**ENT 3.6FL**  
*with Integral Light Guide*



**ENT 3.6FP**  
*with Detachable Light Guide*

Our ENT 3.6F series flexible nasopharyngoscopes offer high-resolution image quality by virtue of their superior optics, illumination and wide field of view, packaged in a lightweight, sleek design that is easy and convenient to use. They are ideal for a variety of minimally invasive diagnostic and surgical applications. The ENT 3.6F is available in two variations to meet your clinical needs.



## Features & Benefits

- Small, 3.6 mm diameter shaft for enhanced patient comfort
- High articulation range for increased maneuverability
- Superior optics and illumination for accurate, consistent imaging
- High-resolution fiberoptic image with a wide field of view
- Lightweight, sleek, ergonomic design
- Fully immersible, leak testable and video compatible

## Specifications

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	3.6mm
Weight	11.6 oz
Articulation Range	135 degrees up / down
Depth of Field	5mm to 50mm
Field of View	70 degrees
Battery	Lithium Ion



64 Technology Park Road, Sturbridge, MA 01566  
Toll-free: 800.225.7486 Worldwide: 1.508.347.5100  
Fax: 508.347.2380  
Email: sales@optim-llc.com  
Website: www.optim-llc.com

Made in U.S.A.



Optim is ISO 9001:2008 / ISO 13485:2003 Certified