



The ENTity VL Endoscopy System

Flexible Laryngoscope | Digital Video Processor | Imaging Capture Suite and Display

Practical Video Endoscopy System for ENT Diagnostics

Maximize voice and swallowing procedures and the patient experience with chip-on-the-tip flexible laryngoscopy

Simplified Work-flow and Patient Management

Easy-to-use patient exam documentation software with image and video capture for review and archive

Compact and Expandable System

Small, versatile footprint with one-touch camera addition for rigid endoscopy and stroboscopy

ENTity VL Video Laryngoscope

- Distal chip video imaging
- Integrated LED illumination and auto white balance
- Durable design and construction

ENTity VL Imaging Capture Suite and Display

- 22" Touchscreen, Medical Grade HD LED Panel
- Anti-microbial, Sealed Front Bezel
- Streamlined processor with synchronous outputs
- Optimized video image or full screen presentation
- One-touch image and video capture and review
- Frame-by-Frame and slow motion review
- Simple administration and saving of diagnostic exams
- One-touch switch from video laryngoscope to add-on camera
- Rapid practice-branded report creation to PDF for print or EMR
- HDMI outputs for additional displays

ENTity VL Video Stroboscopy

- True "flashlight" stroboscopy
- Rolling "Slow Motion" phase and "Fixed" phase options
- Silent operation and compact footprint
- Convenient foot pedal controls



The ENTity VL Endoscopy & Stroboscopy System



ENTity VL Video System

- Video Laryngoscope
- Digital Video Processor
- Imaging Capture Suite and Display
- Endoscope Leak Tester

ENTity VL Stroboscopy Kit

- Stroboscopy Unit
- 70°, 10mm Rigid Laryngoscope
- External LED Light Source for Rigid
- VidCap Add-on Camera System

ENTity VL Peripherals & Accessories

- Transport Cart
- Imaging Capture Suite Stand (non-cart)
- Scope mount for Transport Cart

Includes ENTity System standard 2 Year Warranty

