



The ENTity VL Endoscopy System

Practical Video Endoscopy System for SLP Diagnostics

Maximize voice and swallowing procedures and the patient experience with chip-on-the-tip flexible laryngoscopy

Streamlined Stroboscopy

Rigid and flexible endoscope options with true “flashlight” LED strobing technology

Simplified Work-flow and Patient Management

Easy-to-use patient exam documentation software with image and video capture for review and archive

Compact and Expandable System

Small, versatile footprint with one-touch camera addition for rigid endoscopy and stroboscopy

ENTity VL Video Laryngoscope

- Distal chip video imaging
- Ergonomic design for comfort
- Integrated LED illumination and auto white balance

ENTity VL Imaging Capture Suite and Display

- 22" Touchscreen, anti-microbial medical grade HD LED Panel
- Optimized video image or full screen presentation
- One-touch image and video capture and review
- Frame-by-Frame and slow motion review
- Simple administration and saving of diagnostic exams
- One-touch switch from video laryngoscope to add-on camera
- Rapid practice-branded report creation to PDF for print or EMR

ENTity VL Video Stroboscopy

- True “flashlight” stroboscopy
- Rolling “Slow Motion” phase and “Fixed” phase options
- Silent operation and compact footprint
- Rigid and flexible options
- Convenient foot pedal controls



The ENTity VL Endoscopy System



ENTity VL Video System

- Video Laryngoscope
- Digital Video Processor
- Imaging Capture Suite and Display
- Endoscope Leak Tester

ENTity VL Stroboscopy Kit

- Stroboscopy Unit
- Rigid Laryngoscope
- External LED Light Source for Rigid
- VidCap Add-on Camera System

ENTity VL Peripherals & Accessories

- Transport Cart
- Imaging Capture Suite Stand (non-cart)
- Scope mount for Transport Cart

Includes ENTity System standard 2 Year Warranty

