

Product
EXCLUSIVE

EndoMersion™

Full Endoscope Immersion System

Features & Benefits

- **Fully Immerse Endoscope**— Completely submerge the entire endoscope for HLD processing for patient and staff safety
- **Vertical Wall Mount**— Convenient staff access while vertical position keeps optical insertion shaft in ideal position for protection and storage
- **Unique Design**— Endoscope lifter with lid provides endoscope security for safe handling and minimal staff exposure to disinfectant solutions
- **Ease of Use**— Transparent chambers allow for instant inspection and economical 1½ quart solution volume

Dimensions:

Wall Bracket | 24"L x 5.5"W x 5"D

Tube | 28"L



64 Technology Park Road, Sturbridge, MA 01566 USA
Phone: +1.508.347.5100 | Fax: +1.508.347.2380
sales@optim-llc.com | www.optim-llc.com

