

# Flexible Nasopharyngoscope Model ENT 3.6FP

*Precision endoscopy with proven performance*



## ENT 3.6FP with Detachable Light Guide

Our ENT 3.6FP series Flexible Nasopharyngoscope offers high-resolution image quality by virtue of their superior optics, illumination and wide field of view, packaged in a lightweight, sleek design that is easy and convenient to use. Ideal for a variety of minimally invasive diagnostic and surgical applications. The ENT 3.6FP is available to meet your clinical needs.

The ENT 3.6FP detachable light guide connection gives the physician options for connecting to a light source. The Optim ENT 3.6FP is compatible with the Optim Portable LED Light Source. This portable solution eliminates the need for external light sources and lamp replacement and allows for the freedom to move from room to room.



### Features & Benefits

- 3.6 mm diameter insertion tube for enhanced patient comfort
- High articulation range for increased maneuverability
- Superior optics and illumination for accurate, consistent imaging
- High-resolution fiberoptic image with a wide field of view
- Lightweight, sleek, ergonomic design
- Fully immersible, leak testable and video compatible

### Specifications

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	3.6mm
Weight	11.6 oz
Articulation Range	135 degrees up / down
Depth of Field	5mm to 50mm
Field of View	70 degrees
Battery	Lithium Ion

Contact us for more information on the Portable LED Light Source



64 Technology Park Road, Sturbridge, MA 01566 USA  
Phone: +1.508.347.5100 | Fax: +1.508.347.2380  
sales@optim-llc.com | www.optim-llc.com



Optim is ISO13485:2016 Certified