

The ENTity Series of Flexible Fiberoptic Nasopharyngoscopes can be used virtually anywhere you need to be. It is the only flexible endoscope that carries a guaranteed two (2) year warranty due to its reliability and durability. Our patented OptiLux integrated LED illumination system eliminates cumbersome external light boxes while providing exceptional visibility, making the endoscope completely portable. The OptiLux lighting system within the endoscope handle provides superior lighting targeted in the 5500K color temperature range for superior tissue color differentiation.



**Industry's only Nasopharyngoscope with Self-Contained Light Source & 2-Year Warranty**

**ENTITY**  
SDXL

ENTity SDXL features a slimline 2.7mm insertion shaft making it the preferred choice for diagnostic evaluation of narrow upper airway passages in both adults and children.

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	2.7mm
Weight	11.6 oz
End Tip Articulation	110 degrees up / down
Depth of Field	5mm to 50mm
Field of View	85 degrees
Battery	Lithium-Ion

**ENTITY**  
XL

ENTity XL makes diagnostic examination of the nose and throat for viewing of the nasal passageways and larynx as easy to use as a stethoscope.

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	3.6mm
Weight	11.6 oz
End Tip Articulation	135 degrees up / down
Depth of Field	5mm to 50mm
Field of View	70 degrees
Battery	Lithium-Ion

**Better economics**

- See more patients per day, increase revenue
- Reduced capital expenses and ongoing maintenance
- No lamp or light guide replacement costs

**Portable, compact, all-in-one design**

- No external light source required
- Fully sealed and completely soakable
- Ideal for office, bedside & hospital use