

ENT Endoscope Storage Cabinet

Stainless steel, 6 ENT endoscope capacity, wall mount storage cabinet with key lock



Features & Benefits

- Customizable endoscope holders to fit most common endoscope brands in the field
- Stores up to 6 ENT endoscopes
- Glass paneled doors for easy viewing of endoscope inventory
- Industry best 12 year limited manufacturer warranty
- Ideal for endoscopes without integral light guide cable or removable video/light processor cord
- Cabinet dimensions: 30”H x 24”W x 12”D
- Cabinet weight: 60 lbs

Optional LED-UV disinfection strip helps to keep scopes disinfected!



HEPA Filter Option

HEPA Filtered Positive Pressure with two high volume intake fans and one exhaust fan that activate automatically upon closing of doors. Change HEPA filters without opening cabinet, save time and ensure no cycle interruptions



64 Technology Park Road, Sturbridge, MA 01566 USA
Phone: +1.508.347.5100 | Fax: +1.508.347.2380
sales@optim-llc.com | www.optim-llc.com



Optim is ISO13485:2016 Certified