

Rigid Sinuscope & PLS Bundle

The ENT Rigid Sinuscopes are designed and manufactured to the highest quality standards ensuring reliable performance to meet the rigors of real-world office and surgical use. The rigid sinuscopes combine precision machined components with superior optics to provide exceptional clarity and brightness. When combined with the Portable LED Light Source (PLS) it creates a portable solution eliminating the need for external light sources and lamp replacement.



Specifications

Diameter	Working Length	Viewing Direction	Field of View
4.0 mm	175 mm	0°	standard
4.0 mm	175 mm	30°	standard
3.0 mm	175 mm	0°	standard
3.0 mm	175 mm	30°	standard
2.7 mm	175 mm	0°	standard
2.7 mm	175 mm	30°	standard

Optim ENT Rigid Sinuscopes are soakable, low temp autoclave and Sterrad compatible

Features & Benefits


- **Portability**—ENT Rigid Sinuscope with PLS eliminates cumbersome cords, providing freedom of movement during an exam
- **Compatibility**—Universal eyepiece is compatible with add-on camera systems
- **Power Source**—Rechargeable batteries provide power for an entire day of use
- **Variable Light Intensity**—Adjustable light control intensity switch
- **Integrated Lighting Design**—Our patented light coupling technology, OptiLux focuses and shapes the light while delivering continuous light intensity where it counts: at the end of the scope

See for yourself how the Rigid/PLS Bundle can change your daily procedures



Rigid Sinuscope & PLS Bundle

*Bundle includes:
(1) Rigid Sinuscope, Portable LED Light Source with (2) batteries and charger in a compact carrying case.*



Rigid Sinuscope Model 4.0mm 175mm with Portable Light Source



64 Technology Park Road, Sturbridge, MA 01566 USA
 Phone: +1.508.347.5100 | Fax: +1.508.347.2380
 sales@optim-llc.com | www.optim-llc.com



Optim is ISO13485:2016 Certified