

Endoscope Storage Caddies



Part No. F004732P (EndoCaddy® AZ2020)

Single tube holder for flexible nasopharyngoscopes, complete with one 20" long x 1 1/4" diameter tube, bracket, vapor cap and mounting bolts

Part No. F005534P (EndoCaddy® AZ4040)

Double tube holder for flexible nasopharyngoscopes, complete with two 20" long x 1 1/4" diameter tubes, bracket, vapor caps and mounting bolts

Part No. F014395P (EndoCaddy® AZ4343)

Triple tube holder for flexible nasopharyngoscopes, complete with three 20" long x 1 1/4" diameter tubes, bracket, vapor caps and mounting bolts

Part No. F005533P (EndoCaddy® AZ5050)

Four tube holder for flexible and rigid endoscopes, complete with two 20" long x 1 1/4" and two 14" long x 3/4" tubes, bracket, vapor caps and mounting bolts