

## ENT 3.6FP Flexible Nasopharyngoscope

*Precision endoscopy with proven performance*

Our ENT 3.6FP series Flexible Nasopharyngoscope offers high-resolution image quality by virtue of their superior optics, illumination and wide field of view, packaged in a lightweight, sleek design that is easy and convenient to use. Ideal for a variety of minimally invasive diagnostic applications.

The ENT 3.6FP detachable light guide connection provides options for connecting to a light source. The ENT 3.6FP is compatible with the Optim Portable LED Light Source. This portable solution eliminates the need for external light sources and lamp replacement and allows for the freedom to move from room to room.

**Contact us today!**

+1.508.347.5100

sales@optim-llc.com

### Specifications

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	3.6mm
Weight	11.6 oz
Articulation Range	135 degrees up / down
Depth of Field	5mm to 50mm
Field of View	70 degrees
Light Guide Connection	Detachable

### Kit Includes

ENT 3.6FP Endoscope, Storz & Wolf Light Guide Adapters, Vent Tag & Instruction Card, IFU and Carrying Case



### Features & Benefits

**Patient Comfort**—3.6mm insertion tube diameter for enhanced patient comfort.

**Quality Optics**— Superior resolution fiber optics and illumination for accurate, consistent imaging and diagnosis.

**Ease of Use**—Lightweight, ergonomic handle design with high articulation range for increased maneuverability.

**Optional Lighting**— Post style light guide connection provides options for connecting to a light source.

**Video Compatibility**—Compatible with HD VidCap USB Camera and TIMS Software.



The ENT 3.6FP is compatible with the Optim Portable LED Light Source. This portable solution eliminates the need for external light sources and lamp replacement.



The ENT 3.6FP allows for connection to standard light source box utilizing standard light guide cables enabling connectivity to any light source on the market.



64 Technology Park Road  
Sturbridge, MA 01566 USA

+1.508.347.5100  
sales@optim-llc.com

**OPTIM'S PRODUCTS ARE  
PROUDLY MANUFACTURED  
IN THE USA**