

ENTity MobileFEES | Portable On-The-Go Kit

Leading Performance. Optimal Simplicity. Proven Support

Portable On-The-Go Swallowing Diagnostics
Examine | Evaluate | Diagnose



The MobileFEES (Flexible Endoscopic Evaluation of Swallowing) Portable On-the-Go Kit, allows Clinicians the choice of bedside portability or office based FEES programs. The Optim FEES Kit is the optimal choice for nursing homes, rehabilitation settings, hospitals, clinics and independent SLP practices.

Contact us today!

+1.508.347.5100

sales@optim-llc.com



MobileFEES | Portable On-the-Go Kit:

ENTity Flexible Nasopharyngoscope w/ Internal LED Light Source
The ENTity is available in the standard 3.6mm insertion tube diameter as well the small diameter 2.7mm.

HD VidCap USB Video Camera Kit

TIMS Lite or TIMS Dicom Software Available

Software Support and Maintenance Available

Equipment In-Service with our Trained Clinical Specialists

Two-Year Manufacturers Warranty on ENTity XL / SDXL with Free Loaner Equipment

Kit Accessories Include:

Battery Charger, (4) Rechargeable Lithium Ion Batteries, Leak Tester, Foot Pedal Camera Activator & Portable Carrying Case



Training & Service

- Equipment in-service by trained SLP's
- Superior customer service

Portability

- Lightweight, portable case
- Bring the scope to the patient
- Use in the hospital / bedside / office

Simple Software

- Capture still images & streaming videos
- EMR compatible & HIPAA compliant
- Windows based

Industry's only Nasopharyngoscope with Self Contained Light Source & 2-Year Warranty



64 Technology Park Road
Sturbridge, MA 01566 USA

+1.508.347.5100
sales@optim-llc.com

OPTIM'S PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA