

ENTity Flexible Fiberoptic Nasopharyngoscope

Designed to empower physicians with the visibility & mobility they need



Industry's only Nasopharyngoscope with Self-Contained Light Source & 2-Year Warranty

The ENTity Series of Flexible Fiberoptic Nasopharyngoscopes revolutionizes the field of nasal visualization with its cutting-edge features and unparalleled portability. Equipped with the patented OptiLux integrated LED illumination system, the ENTity eliminates the need for bulky external light boxes and cables, allowing for seamless mobility and ease of use in any setting.

Contact us today!

+1.508.347.5100

sales@optim-llc.com

ENTITY SDXL

ENTity SDXL features a slimline 2.7mm insertion shaft making it the preferred choice for diagnostic evaluation of narrow upper airway passages.

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	2.7mm
Weight	11.6 oz
End Tip Articulation	110 degrees up / down
Depth of Field	5mm to 50mm
Field of View	85 degrees
Battery	Lithium-Ion

Kit Includes

ENTity Endoscope, (2) Li-Ion Batteries, Battery Charger, Vent Tag & Instruction Card, IFU and Carrying Case



ENTITY XL

ENTity XL makes diagnostic examination of the nose and throat for viewing of the nasal passageways and larynx as easy to use as a stethoscope.

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	3.6mm
Weight	11.6 oz
End Tip Articulation	135 degrees up / down
Depth of Field	5mm to 50mm
Field of View	70 degrees
Battery	Lithium-Ion

Features & Benefits

Better Economics—See more patients per day, increasing revenue. Reduce capital expenses and ongoing maintenance - no lamp or light guide replacement costs.

Portability—Bring the scope to the patient; ideal for office, bedside & hospital use.

Patented Lighting—Powered by OptiLux Lighting, ENTity provides superior lighting targeted at 5500K color temperature for superior tissue color differentiation.

Video Compatibility—Compatible with HD VidCap USB Camera and TIMS Software.



64 Technology Park Road
Sturbridge, MA 01566 USA

+1.508.347.5100
sales@optim-llc.com

OPTIM'S PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA