

The system's upright wall mount design allows for easy access by staff, promoting streamlined workflow and quick deployment when needed. Its vertical orientation ensures that the insertion tube is kept in the ideal position for protection and storage, minimizing the risk of damage or deterioration.

EndoCaddy® AZ4040, double tube holder is ideal for soaking an endoscope while utilizing the other tube for storage or drying of another flexible endoscopes, Nasopharyngoscope or Rhino-Laryngoscope.

Easily mounts to wall, cabinet or cart to protect or allow soaking of ENT flexible endoscope. Storage tube can be removed from bracket for dispensing solutions or cleaning.

Contact us today!

+1.508.347.5100

sales@optim-llc.com



**Part No. F004732P (EndoCaddy® AZ2020)**

Single tube holder for flexible nasopharyngoscopes, complete with one 20" long x 1 1/4" diameter tube, bracket, vapor cap and mounting bolts



**Part No. F005534P (EndoCaddy® AZ4040)**

Double tube holder for flexible nasopharyngoscopes, complete with two 20" long x 1 1/4" diameter tubes, bracket, vapor caps and mounting bolts



**Part No. F014395P (EndoCaddy® AZ4343)**

Triple tube holder for flexible nasopharyngoscopes, complete with three 20" long x 1 1/4" diameter tubes, bracket, vapor caps and mounting bolts



**Part No. F005533P (EndoCaddy® AZ5050)**

Four tube holder for (2) flexible and (2) rigid endoscopes, complete with two 20" long x 1 1/4" and two 14" long x 3/4" tubes, bracket, vapor caps and mounting bolts

