

HD VidCap USB Camera with TIMS Patient Management Software

Designed to empower physicians with the visibility & mobility they need

The add-on HD VidCap USB Camera System makes recording video and taking still images easy and convenient. Simply couple to your flexible or rigid endoscope, plug the camera into the USB port on your PC or laptop and you are ready to go. It adds a new dimension of functionality to your endoscope, aiding in diagnosis by colleagues.

The easy to use, reliable and efficient TIMS Patient Management software includes patient database for visit comparison, image annotation, report creation for easy sharing and audio and video recording.



Contact us today!

+1.508.347.5100

sales@optim-llc.com

TIMS Patient Management Software

Compatibility—Record studies at 30fps in full color at maximum resolution with synced audio

Study Storage Options—Archive studies to any PACS, VNA and EMR with the FULL TIMS Version or store directly to device with the TIMS LITE Version.

Productivity— Quickly review video/image for annotation, reporting and more

Intuitive Software—Eliminate the out of date workflows and software. TIMS is simple to use and includes powerful editing and analysis tools

HIPPA Compliant—Keep workflow digital. No memory sticks, CDs or DVDs to store and transport files



Optim's ENTity XL Nasopharyngoscope coupled to the HD VidCap Camera then connected to a laptop for easy image viewing.

HD VidCap USB Plug and Play Camera

Easy to use—USB Plug & Play, intuitive design, compatible with any qualified Windows-based computer, automatic anti-moiré filter and white balance

Compatibility—Compatible with flexible and rigid endoscopes

Patient and Office Safety—Allows for separation between you and the patient during exam

Disinfection—Camera and coupler are fully immersible

Sensor	CMOS 1/2" high sensitivity
Resolution	1920 x 1080
Weight	210g
Cable Length	3 meters
Protection Type	IPX7 (waterproof)
Additional Specs	Multi function buttons for photo and movie capture, digital zoom, white balance, anti-moiré filter

PC/Laptop Minimum Specifications for TIMS Software:

- 10th Generation Intel Core i5 Processor
- 16GB Memory
- 500GB SSD
- Windows 10 or 11 Pro



64 Technology Park Road
Sturbridge, MA 01566 USA

+1.508.347.5100
sales@optim-llc.com

OPTIM'S PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA