


Rigid Sinusscopes & PLS Bundle

High Quality, With A Versatile Portable Design For Ease Of Use



The ENT Rigid Sinusscopes are designed and manufactured to the highest quality standards ensuring reliable performance to meet the rigors of real-world office and surgical use. The rigid sinusscopes combine precision machined components with superior optics to provide exceptional clarity and brightness. When combined with the Portable LED Light Source (PLS) it creates a portable solution eliminating the need for external light sources and lamp replacement.

Specifications

| Diameter | Working Length | Viewing Direction | Field of View |
|----------|----------------|-------------------|---------------|
| 4.0 mm | 175 mm | 0° | standard |
| 4.0 mm | 175 mm | 30° | standard |
| 3.0 mm | 175 mm | 0° | standard |
| 3.0 mm | 175 mm | 30° | standard |
| 2.7 mm | 175 mm | 0° | standard |
| 2.7 mm | 175 mm | 30° | standard |

Contact us today!

+1.508.347.5100

sales@optim-llc.com

See for yourself how the Rigid/PLS Bundle can change your daily procedures



Features & Benefits

Portability—ENT Rigid Sinuscope with PLS eliminates cumbersome cords, providing freedom of movement during an exam

Compatibility—Universal eyepiece is compatible with add-on camera systems

HLD—Rigid Sinusscopes are soakable, low temp autoclave and Sterrad compatible

Power Source—Rechargeable batteries provide power for an entire day of use

Kit Includes

(1) Rigid Sinuscope, Portable LED Light Source with (2) batteries and charger in a compact carrying case.



64 Technology Park Road
Sturbridge, MA 01566 USA

+1.508.347.5100

sales@optim-llc.com

OPTIM'S PRODUCTS ARE
PROUDLY MANUFACTURED
IN THE USA